

DMEM高糖培养基相关产品成分对比表								
货号	L-谷氨酰胺	L-丙氨酰-L-谷氨酰胺	丙酮酸钠	HEPES	酚红	双抗	Ca ²⁺	Mg ²⁺
EC11078	√		√		√		√	√
EC11078A	√		√		√	√	√	√
EC11069							√	√
EC11077	√				√		√	√
EC11077A	√				√	√	√	√
EC11295			√		√		√	√
EC11295A			√		√	√	√	√
EC11296					√		√	√
EC11116	√			√	√		√	√
EC11116A	√			√	√	√	√	√
EC11115	√		√	√	√		√	√
EC11297			√	√	√		√	√
EC11298				√	√		√	√
EC11299	√		√				√	√
EC11299A	√		√			√	√	√
EC11300			√				√	√
EC11301	√						√	√
EC11302	√		√	√			√	√
EC11303			√	√			√	√
EC11304	√			√			√	√
EC11305				√			√	√
EC11306		√	√		√		√	√
EC11307		√	√				√	√
EC30022		√			√		√	√
EC11308		√					√	√
EC11309		√	√	√	√		√	√
EC11310		√		√	√		√	√
EC11311		√	√	√			√	√
EC11312		√		√			√	√

北京伊事达科技有限公司
电话:13564444959
官网:www.followme-shop.com
地址:北京市海淀区东北旺西路 58 号尚科办公社区 C 区一楼



公众号



客服